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

Oxygen
Concentration
Of Media

Aerobic and anaerobic **bacteria** can be identified by growing them in test tubes of **thioglycollate broth**:

- 1: **Obligate aerobes** need oxygen because they cannot ferment or respire anaerobically. They gather at the top of the tube where the oxygen concentration is highest.
- 2: **Obligate anaerobes** are poisoned by oxygen, so they gather at the bottom of the tube where the oxygen concentration is lowest.
- 3: **Facultative anaerobes** can grow with or without oxygen because they can metabolise energy aerobically or anaerobically. They gather mostly at the top because aerobic respiration generates more ATP than either fermentation or anaerobic respiration.

- 4: **Microaerophiles** need oxygen because they cannot ferment or respire anaerobically. However, they are poisoned by high concentrations of oxygen. They gather in the upper part of the test tube but not the very top.
- 5: **Aerotolerant organisms** do not require oxygen as they metabolise energy anaerobically. Unlike obligate anaerobes however, they are not poisoned by **oxygen**. They can be found evenly spread throughout the test tube.

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1.1 재료 : **Standard Methods Agar (Plate Count Agar)**

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Plate Count Agar (A.P.H.A.)

(Standard Methods Agar, Tryptone Glucose Yeast Agar)

NCM0010A / 500g

Formula	Liter
Tryptone	5.0g
Yeast Extract	2.5g
Glucose	1.0g
Agar No.1	15.0g

Description

Formulated to A.P.H.A. specifications, this medium is used for establishing total viable counts for aerobes in food, dairy and water bacteriology. The product uses agar of very high gel strength in order that it can be used in pour-plate as well as surface inoculation techniques. The product can be remelted prior to use although it should not be held for a prolonged period in the molten stage.

Method for reconstitution

Weigh 23.5 grams of powder, disperse in 1 litre of deionised water. Bring to the boil with frequent stirring to dissolve. Dispense into tubes and sterilise by autoclaving at 121°C for 15 minutes. Cool to 44-46°C for not more than 3 hours prior to use. ROLL-TUBES. Add an additional 10g/litre Agar No. 1 prior to reconstitution of the medium.

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